


ARG83750 Human Perforin ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83750 Human Perforin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Perforin in serum, plasma and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Perforin
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	46.88 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	93.75 - 6000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83750 Human Perforin ELISA Kit standard curve image

ARG83750 Human Perforin ELISA Kit results of a typical standard run with optical density reading at 450 nm.