


## ARG83755 Human RANTES ELISA Kit

Package: 96 wells  
Store at: 4°C, -20°C

### Summary

|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Description | ARG83755 Human RANTES ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human RANTES in serum, plasma and cell culture supernatants. |
| Tested Reactivity   | Hu                                                                                                                                                  |
| Tested Application  | ELISA                                                                                                                                               |
| Target Name         | RANTES                                                                                                                                              |
| Conjugation         | HRP                                                                                                                                                 |
| Conjugation Note    | Substrate: TMB and read at 450 nm.                                                                                                                  |
| Sensitivity         | 7.81 pg/ml                                                                                                                                          |
| Sample Type         | Serum, plasma and cell culture supernatants.                                                                                                        |
| Standard Range      | 15.63 - 1000 pg/ml                                                                                                                                  |
| Sample Volume       | 50 µl                                                                                                                                               |
| Precision           |                                                                                                                                                     |

### Images



ARG83755 Human RANTES ELISA Kit standard curve image

ARG83755 Human RANTES ELISA Kit results of a typical standard run with optical density reading at 450 nm.