


ARG83758 Mouse FGF basic ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83758 Mouse FGF basic ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse FGF basic in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	FGF basic
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	23.44 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	46.88 - 3000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83758 Mouse FGF basic ELISA Kit standard curve image

ARG83758 Mouse FGF basic ELISA Kit results of a typical standard run with optical density reading at 450 nm.