


ARG83759 Mouse CXCL1 / GRO alpha ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83759 Mouse CXCL1 / GRO alpha ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse CXCL1 / GRO alpha in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	CXCL1 / GRO alpha
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	7.81 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	15.63 - 1000 pg/ml
Sample Volume	50 µl
Precision	

Images



ARG83759 Mouse CXCL1 / GRO alpha ELISA Kit standard curve image

ARG83759 Mouse CXCL1 / GRO alpha ELISA Kit results of a typical standard run with optical density reading at 450 nm.