

## ARG83808 beta-Glucosidase Activity Assay Kit (Colorimetric)

Package: 100 tests

Store at: -20°C

### Summary

Product Description	ARG83808 beta-Glucosidase Activity Assay Kit (Colorimetric) can be used to measure beta-Glucosidase activity in biological samples.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	beta-Glucosidase
Conjugation Note	Read at 405 nm
Sensitivity	2 U/L
Detection Range	2 - 250 U/L
Sample Type	biological samples.
Sample Volume	20 µl

### Application Instructions


Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

### Properties

Form	Liquid
Storage instruction	Store components at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Symbol	GBA
Gene Full Name	beta-Glucosidase



ARG83808 beta Glucosidase Activity Assay Kit (Colorimetric) enzyme kinetics graph

Kinetics of 0 and 250 U/L beta Glucosidase reaction, using ARG83808 beta Glucosidase Activity Assay Kit (Colorimetric).