

ARG83813
DDDDK tag ELISA KitPackage: 96 wells
Store at: 4°C**Summary**


Product Description	ARG83813 DDDDK tag ELISA Kit is an Enzyme Immunoassay kit for the quantification of DDDDK tag in serum, plasma and cell culture supernatants.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	DDDDK tag
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	15.6 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	31.3-1000 ng/ml
Sample Volume	100 µL

Application Instructions

Assay Time	2 hours
------------	---------

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG83813 DDDDK tag ELISA Kit standard curve image

ARG83813 DDDDK tag ELISA Kit results of a typical standard run with optical density reading at 450 nm.