

ARG83838 Glutamate Assay Kit (Colorimetric)

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83838 Glutamate Assay Kit (Colorimetric) can be used to measure Glutamate in urine, serum, plasma, tissue extracts, cell lysate, cell culture media, other biological fluids.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Glutamate
Conjugation Note	Read at 450 nm
Sample Type	urine, serum, plasma, tissue extracts, cell lysate, cell culture media, other biological fluids
Standard Range	0.03 mmol/L - 3 mmol/L
Precision	10 µl

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.