

ARG83839
Phenylalanine Ammonia-Lyase Activity Assay Kit (Colorimetric)

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG83839 Phenylalanine Ammonia-Lyase Activity Assay Kit (Colorimetric) can be used to measure Phenylalanine Ammonia-Lyase Activity in tissue extracts, cell lysate.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Phenylalanine Ammonia-Lyase
Conjugation Note	Read at 620 nm
Sample Type	tissue extracts, cell lysate
Standard Range	10 µmol/L - 1000 µmol/L
Precision	10 µl

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.