

ARG83840
Na⁺ / K⁺ ATPase Activity Assay Kit (Colorimetric)

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG83840 Na ⁺ / K ⁺ ATPase Activity Assay Kit (Colorimetric) can be used to measure Na ⁺ / K ⁺ ATPase Activity in tissue extracts and cell lysate.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Na ⁺ / K ⁺ ATPase
Conjugation Note	Read at 650 nm
Sample Type	Tissue extracts, Cell lysate
Standard Range	10 µmol/L - 1000 µmol/L
Sample Volume	10 µl

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.