

ARG83853
PDHX / Pyruvate dehydrogenase Activity Assay Kit (Colorimetric)

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG83853 PDHX / Pyruvate dehydrogenase Activity Assay Kit (Colorimetric) can be used to measure PDHX / Pyruvate dehydrogenase Activity in tissue, cultured cell, and mitochondria isolates.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	PDHX / Pyruvate dehydrogenase
Conjugation Note	Read at 450 nm
Sample Type	Tissue, Cultured Cell, Mitochondria isolates
Standard Range	12.5 µmol/L - 400 µmol/L
Sample Volume	20 µl

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.