

ARG83876
Cell Counting Kit-8 Assay KitPackage: 100 tests
Store at: 4°C

Summary

Product Description	ARG83876 Cell Counting Kit-8 Assay Kit is used for the quantitation of viable cell numbers in proliferation and cytotoxicity assays.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Cell Counting Kit-8
Conjugation Note	Read at 450 nm

Application Instructions

Application Note	Please note that this kit does not include a microplate.
------------------	--

Properties

Form	Liquid
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.