

ARG83877 Mitochondrial Membrane Potential (JC-1) Assay Kit

Package: 100 tests

Store at: -20°C

Summary

Product Description	ARG83877 Mitochondrial Membrane Potential (JC-1) Assay Kit is used for the detection of mitochondrial membrane potential using a fluorescence microplate reader.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Mitochondrial Membrane Potential (JC-1)
Conjugation Note	λ ex/em = 515/529; 585/590 nm (fluorometric)

Application Instructions

Application Note	Please note that this kit does not include a microplate.
------------------	--

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.