

ARG83878
Intracellular ROS Assay Kit (Fluorometric)Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83878 Intracellular ROS Assay Kit (Fluorometric) can be used to measure intracellular Reactive Oxygen Species (ROS) levels in suspension or adherent cell cultures.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Intracellular ROS
Conjugation Note	λ ex/em = 495/529 nm (fluorometric)
Sample Type	suspension or adherent cell cultures

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.