

ARG83880 Human Vasorin ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83880 Human Vasorin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Vasorin in tissue homogenates, cell lysates, biological fluids
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Vasorin
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.16 ng/mL
Sample Type	tissue homogenates, cell lysates, biological fluids
Standard Range	0.32-20 ng/mL
Sample Volume	100 µL
Alternate Names	SLITL2; Virga; SLIT-Like 2

Application Instructions


Assay Time	3.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	VASN
Gene Full Name	Vasorin
Background	Enables transforming growth factor beta binding activity. Involved in negative regulation of epithelial to mesenchymal transition and negative regulation of transforming growth factor beta receptor signaling pathway. Located in cell surface and extracellular exosome. [provided by Alliance of Genome Resources, Jul 2025]
Function	Acts as an inhibitor of TGF-beta signaling. [UniProt]
PTM	Disulfide bond; Glycoprotein. [UniProt]
Cellular Localization	Membrane; Secreted. [UniProt]



ARG83880 Human Vasorin ELISA Kit standard curve image

ARG83880 Human Vasorin ELISA Kit results of a typical standard run with optical density reading at 450 nm.