

ARG83882

arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha)

Package: 96 wells

Store at: 4°C, -20°C

Summary

Product Description

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) is an Enzyme Immunoassay kit for the quantification of Human CAR-T associated Cytokine in serum, plasma and cell culture supernatants.

[See all Multiplex ELISA kits](#)

Tested Reactivity

Hu

Tested Application


ELISA

Key features

Images


	1	2	3	4	5	6	7	8	9	10	11	12
A	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ	IFN-γ
B	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β	IL-1β
C	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6	IL-6
D	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8	IL-8
E	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10	IL-10
F	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15	IL-15
G	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF	GM-CSF
H	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α	TNF-α

Antibodies Coating Pattern In Microtiter Plate of ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha)




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IFN gamma run with optical density reading at 450 nm.




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IL1 beta run with optical density reading at 450 nm.




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IL6 run with optical density reading at 450 nm.




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IL8 run with optical density reading at 450 nm.




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IL10 run with optical density reading at 450 nm.




ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human IL15 run with optical density reading at 450 nm.



ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human GM-CSF run with optical density reading at 450 nm.



ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) standard curve image

ARG83882 arigoPLEX® Human CAR-T associated Cytokine Multiplex ELISA Kit (IFN gamma, IL1 beta, IL6, IL8, IL10, IL15, GM-CSF, TNF alpha) results of a typical standard for Human TNF alpha run with optical density reading at 450 nm.