

ARG83883 Mouse TSH ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83883 Mouse TSH ELISA Kit is an enzyme immunoassay kit for the quantification of Mouse TSH (thyrotropin) in serum and plasma (EDTA).
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	TSH
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.3125 ng/ml
Sample Type	Serum and plasma (EDTA).
Standard Range	0.625-10 ng/ml
Sample Volume	10 µl
Alternate Names	TSH-B; TSH-beta; Thyroid-stimulating hormone subunit beta; Thyrotropin subunit beta; Thyrotropin beta chain; Thyrotropin alfa; TSH-BETA

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	TSHB
Gene Full Name	thyroid stimulating hormone, beta
Background	The four human glycoprotein hormones chorionic gonadotropin (CG), luteinizing hormone (LH), follicle stimulating hormone (FSH), and thyroid stimulating hormone (TSH) are dimers consisting of alpha and beta subunits that are associated noncovalently. The alpha subunits of these hormones are identical, however, their beta chains are unique and confer biological specificity. Thyroid stimulating hormone functions in the control of thyroid structure and metabolism. The protein encoded by this gene is the beta subunit of thyroid stimulating hormone. Mutations in this gene are associated with congenital central and secondary hypothyroidism and Hashimoto's thyroiditis. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, May 2013]
Function	Indispensable for the control of thyroid structure and metabolism. [UniProt]
Highlight	Related products: TSH antibodies ; TSH ELISA Kits ; TSH Duos / Panels ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData

