




Autophagy


Antibody • Duo • ELISA kit

Autophagy is a process of self-degradation of cellular components. It is activated in response to cellular stress including starvation, growth factor deprivation, ER stress and pathogen infection. Autophagy also plays a critical role in development and differentiation. It is important not only for balancing energy sources in development or in response to cellular stress but also for removing misfolded/aggregated proteins, damaged organelles and pathogens. Autophagy is associated with human pathologies including cancer, neurodegeneration and infectious diseases.




is better than one
arigo speciality


LC3A and LC3B Antibody Duo (ARG30273)




LC3A antibody (ARG51300)



LC3B antibody (ARG55799)



UVRAG antibody (ARG54897)



ATG4D antibody (ARG54820)

arigo Cat.	Product	Clonality	Reactivity	Application
ARG54816	ATG3 antibody	Mouse mAb	Hu, Ms	WB, IHC-P
ARG54904	ATG4A antibody	Rabbit pAb	Hu, Ms	WB, IP
ARG54818	ATG4B antibody	Rabbit pAb	Hu	ICC/IF, WB, IHC-P
ARG54819	ATG4C antibody	Rabbit pAb	Hu, Ms	WB, IHC-P
ARG54820	ATG4D antibody	Mouse mAb	Hu	WB, ICC/IF
ARG54822	ATG5 antibody	Rabbit pAb	Hu, Ms	ICC/IF, IHC, WB
ARG54821	ATG5 antibody	Mouse mAb	Hu, Ms, Rat	WB
ARG54823	ATG7 antibody	Rabbit pAb	Hu, Ms	ICC/IF, WB, IHC-P
ARG54825	ATG9A antibody	Rabbit pAb	Hu	ICC/IF, WB, IHC-P
ARG54815	ATG12 antibody	Rabbit pAb	Hu, Ms	ICC/IF, WB
ARG55112	ATG13 antibody	Rabbit pAb	Hu, Ms, Rat	ELISA, ICC/IF, IHC, WB
ARG54957	ATG16L antibody	Rabbit pAb	Hu, Ms	ICC/IF, WB, IHC-P
ARG54919	ATG101 antibody	Rabbit pAb	Hu, Ms, Rat	ELISA, IHC, WB, ICC/IF
ARG65614	Beclin-1 antibody	Rabbit pAb	Hu, Ms, Rat	ELISA, IHC-P, WB
ARG54094	Beclin-1 antibody	Mouse mAb	Hu, Ms, Rat, Mk, Hm	WB
ARG54959	Beclin-1 antibody	Mouse mAb	Hu	ICC/IF, WB, IHC-P
ARG52234	Beclin-1 pS234 antibody	Rabbit pAb	Hu	WB
ARG52235	Beclin-1 pS295 antibody	Rabbit pAb	Hu	WB
ARG54833	BNIP3 antibody	Rabbit pAb	Hu, Ms	ICC/IF, WB, IHC-P
ARG55622	LC3 antibody	Mouse mAb	Hu, Ms	ICC/IF, WB, IHC-P
ARG51300	LC3A antibody	Rabbit pAb	Hu, Ms, Rat	ICC/IF, WB
ARG55799	LC3B antibody	Rabbit pAb	Hu, Ms, Rat	ICC/IF, IHC-P, WB
ARG54708	LC3C antibody	Rabbit pAb	Hu	WB, IHC-P
ARG65701	NBR1 antibody	Rabbit pAb	Hu, Ms	WB, IHC-P, IP
ARG54876	PI3 Kinase Class III antibody	Rabbit pAb	Hu	WB, IHC-P, ICC/IF
ARG55040	SQSTM1 / p62 antibody	Rabbit pAb	Hu, Ms	WB
ARG52384	SQSTM1 / p62 pS28 antibody	Rabbit pAb	Hu	WB
ARG52385	SQSTM1 / p62 pT269 / pS272 antibody	Rabbit pAb	Hu	WB
ARG54811	ULK1 antibody	Rabbit pAb	Hu	ICC/IF, WB
ARG54990	ULK3 antibody	Mouse mAb	Hu, Ms, Rat	WB
ARG54949	ULK3 antibody	Rabbit pAb	Hu	ELISA, IHC, WB, ICC/IF
ARG54897	UVRAG antibody	Rabbit pAb	Hu, Ms	ICC/IF, WB, IHC-P



www.arigobio.com
E : info@arigobio.com